

DT-6410 Platinum High-Accuracy Thermometer

The DT-6410 Thermometer displays minimum/maximum temperatures for previous 24 hours (for each hour) and previous 7 days (for each day); or previous 30 days (for each day). Time/date key shows the exact time and date for all minimum and maximum readings.

Range: -100.00 to 199.99° C (-148 to 199.99° F)*
Resolution: 0.01° Accuracy: ±0.1° C

Thermometer is ideal for therapy departments. Triple display shows current ° F/ ° C and minimum/maximum temperatures. Alarm is settable in 0.1° increments and signals when temperature rises above/falls below set points. Visual LEDs and audible alarm signal continuously even if temperature returns to non-alarm range. A 4-wire platinum sensor insures linear response and high accuracy across the entire range. Handle probe model features stainless-steel probe, diameter 0.13 inches, stem length 9.0 inches, overall length 12 inches. Low battery life indicator displays when it's time to change your batteries. The DT-6410 permits easy one-key calibration by user.



Specifications

Range	-100° C to 199.99° C
Accuracy	Within ±0.1° C
Resolution	0.01°
Cable length	3 meters
Display update	1 second
DT-6410 Probe	Stainless steel stem, 3.3 mm diameter x 22.8 cm
Power	Two AAA batteries (included)
Battery life	>1000 hrs. continuous operation
Weight	142 g (5 oz)

* Fahrenheit measurement limited by display