



Table Top Caliper

The table top caliper with a 55 cm scale can only be used to measure from the table top to the top surface of the patient.

Model # 705-660, 55 cm



Body Caliper with Flat Blade

The lightweight body caliper with flat blade is used for measuring body parts. The rigid leg that slips under the patient is 1/8 in thick and 11 in long. The scale is available in 45 cm and 55 cm lengths.

Model # 705-545, 45 cm

Model # 705-550, 55 cm



Under/Over Table Body Caliper

The under/over table body caliper with a 55 cm scale measures 6.5 in from the base plate to the curved blade. The blade is 15.5 in long.

Model # 705-700, 55 cm



Body Caliper with Curved Blade


This lightweight precision body caliper is used for measuring body parts. The rigid leg that slips under the patient is 1.8 in thick and 11 in long. The scale is available in 45 cm and 55 cm lengths.

Model # 705-546, 45 cm

Model # 705-551, 55 cm

Body Calipers

The body calipers are designed with a unique measurement safety feature. The measurement safety is accomplished with one scale and a non-removable slide that prevents errors caused by reading the wrong scale on a dual scale of combination calipers. The slider can be locked at any patient thickness.



Head-Neck Caliper

The head-neck caliper is useful in measurements of either patient diameters or surface to tumor depths. The depth of a tumor or intraoral structure relative to either the lateral or anterior/posterior surfaces can be easily determined. Vertical adjustment is 23 cm, horizontal adjustment is 26 cm.

Model # 705-000, 23 x 26 cm