




Best™ GammaBeam 500 Total Body Irradiator System

The Total Body Irradiator GammaBeam 500 is a teletherapy unit designed to produce a large fixed rectangular radiation field at an extended source-to-skin distance in order to deliver total body irradiation. The unit can also be used for research and dosimetry purposes. Features include: dose rate up to 40 cGy/min at 2.5 m in air, record & verify and imaging capabilities, graphical touch screen data entry interface and motorized vertical motion of the head.



Huestis Medical CVO-2000 Warming Oven for Thermoplastics

- Easy-sliding drawer opens fully to provide maximum access to the material
- Multiple masks can be softened at the same time
- No water — reduces the potential for cross-contamination
- Can be used in CT rooms — no humidity that could otherwise affect the electronics of the imaging system
- Digital temperature setting in Celsius or Fahrenheit
- Insulated case
- Manufactured and assembled in the U.S.A.



Digital dry heat temperature control



You don't go to a shoe store and cut your feet to fit the shoes in stock...then why compromise when it comes to improving the clinical outcome for your patients?

Best Medical is the only company that makes custom seeds and strands to your exact specifications — shipped within 24 hours, 7 days a week, sterile and non-sterile!

Best Iodine-125 Seed



Best Palladium-103 Seed



Best™ Seeds have the best visualization in ultrasound, fluoro, x-ray, CT and MRI, facilitating real-time dosimetry



Best Fiducial Markers