


# Best<sup>TM</sup> medical international

A TEAMBEST GLOBAL COMPANY



**GammaBeam**  
saving lives

## Best GammaBeam<sup>TM</sup> 100/300 Equinox<sup>TM</sup> Teletherapy System with Avanza<sup>TM</sup> 6D Patient Positioning Table

With **NEW** Multi-Leaf Collimator for 80 and 100 cm SAD units—  
IMRT, IGRT, SRS, SBRT and Tomotherapy capable with ActiveRx

### UPGRADE KIT for all old Theratron units, 80 or 100 cm including IMRT capabilities w/built in or external MLC



#### Upgrade includes:

- Removing all old controls, electronics and installing a new control system and covers
- Replacing the old collimator system with the new Equinox collimator
- Replacing the old treatment table with the new Avanza<sup>TM</sup> Table
- Retaining the head rotation capability is optional

#### Upgrade features:

- Calculated Arc Speed
- Graphical Control System
- Asymmetric Jaws (optional)
- R&V System Ready (optional)
- Service Log Files
- On-Board Verification
- Motorized Wedge (optional)
- Collision Detection (optional)



Best<sup>TM</sup> Gammacell<sup>®</sup> 1000/3000  
Blood & Research Irradiator



Best<sup>TM</sup> Raycell  
Mk2 Blood  
& Research  
Irradiator



Best<sup>TM</sup> Raycell Mk1  
Blood & Research  
Irradiator

Best Theratronics Ltd., 413 March Road, Ottawa, Ontario K2K 0E4 Canada  
tel: 613 591 2100 866 792 8598 [www.theratronics.com](http://www.theratronics.com) [www.teambest.com](http://www.teambest.com)

\*Certain products shown are not available for sale currently. TeamBest Global Companies ©2025