

# Best™ nomos® Cyber Sonalis® Ultrasound System



The **Best™ Cyber Sonalis® Ultrasound Imaging System** provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented **SimulView™ Technology** provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.

## Additional Features

- ✓ Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- ✓ 21.5" touch screen technology
- ✓ Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- ✓ Patented probe design
- ✓ Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- ✓ Advanced modular software design provides for future upgrade path
- ✓ Convertible to desktop or portable without losing any functionality
- ✓ Dual battery backup built-in to protect from power loss



Best Nomos • 1 Best Drive • Pittsburgh, PA 15202 USA

tel: 412 312 6700 800 70 NOMOS fax: 412 312 6701 www.nomos.com

\*Certain products shown are not available for sale currently.

# Sonalis® Transducers

We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the **BEST** name in ultrasounds!



## 8L2A Linear Array

APPLICATIONS: Arterial, Carotid, Vascular Access, Venous



## 12L5A Linear Array

APPLICATIONS: Arterial, Breast, Carotid, Dialysis Access, Lung, Neonatal Hip, Nerve Block, Ophthalmic, Testes, Thyroid, Vascular Access, Venous



## 14L3 Linear Array

APPLICATIONS: Arterial, Breast, Carotid, Dialysis Access, Lung, MSK, Neonatal Hip, Nerve Block, Ophthalmic, Testes, Thyroid, Vascular Access, Venous



## 15LW4 Linear Array

APPLICATIONS: Arterial, Breast, Carotid, Dialysis Access, Lung, MSK, Neonatal Hip, Nerve Block, Ophthalmic, Testes, Thyroid, Vascular Access, Venous  
VET BIOPSY KIT AVAILABLE



## 15LA Linear Array

APPLICATIONS: Arterial, Breast, Carotid, Dialysis Access, Lung, MSK, Neonatal Hip, Nerve Block, Ophthalmic, Testes, Thyroid, Vascular Access, Venous  
VET BIOPSY KIT AVAILABLE



## 15L4A Linear Array

APPLICATIONS: Arterial, Breast, Carotid, Dialysis Access, Lung, MSK, Neonatal Hip, Nerve Block, Ophthalmic, Thyroid, Vascular Access, Venous




## 16L5 Linear Array

APPLICATIONS: Breast, Lung, MSK, Nerve Block, Vascular Access  
VET BIOPSY KIT AVAILABLE



## 8V3 Phased Array

APPLICATIONS: Cardiac



## 4V2A Phased Array

APPLICATIONS: Cardiac, FAST, TCD



## 5C2A Curved Array

APPLICATIONS: Abdominal, FAST, Fetal Cardiac, MSK, OB/GYN, Renal, Thyroid, Visceral  
VET BIOPSY KIT AVAILABLE



## 9MC3 Curved Array

APPLICATIONS: Abdominal, Cardiac, Neonatal Head, Small Parts, Thyroid, Vascular Access



## 8EC4A Endocavity

APPLICATIONS: OB/GYN, Prostate  
VET BIOPSY KIT AVAILABLE



## 10EC4A Endocavity

APPLICATIONS: OB/GYN, Prostate  
VET BIOPSY KIT AVAILABLE



## XY-BI-Plane Phased Array

APPLICATIONS: Cardiac, Vascular, Lung



## 10BP4 Bi-Plane

APPLICATIONS: Prostate



## 8BP4 Bi-Plane

APPLICATIONS: Prostate



## 8TE3 Trans-esophageal

APPLICATIONS: Motorized Adult Multiplane TEE Probe



## Pedoff

APPLICATIONS: Cardiac



## 16HL7 High Frequency Linear Array

APPLICATIONS: MSK, Venous



[www.globalbestcure.org](http://www.globalbestcure.org)



[www.bestcure.md](http://www.bestcure.md)



[www.bestcure.md](http://www.bestcure.md)



[www.brachytherapy.org](http://www.brachytherapy.org)

Best Medical International • 7643 Fullerton Road • Springfield, VA 22153  
tel: 703 451 2378 800 336 4970 [www.besttotalsolutions.com](http://www.besttotalsolutions.com) [www.teambest.com](http://www.teambest.com)