


# Best<sup>®</sup> medical international

## Best<sup>®</sup> Radioactive Seeds for Brachytherapy & Breast Cancer Localization


- We manufacture both I-125 and Pd-103 seeds in the U.S. and supply directly to hospitals
- Best Medical is the only company that makes custom seeds and strands to your exact specifications — any source, any configuration, any time!
- Best<sup>®</sup> seeds have the best visualization in ultrasound, fluoro, X-ray, CT and MRI, facilitating real-time dosimetry
- We help set up new programs with licensing, procurement of necessary equipment, proctoring, reimbursement consulting, patient education and marketing



Best<sup>®</sup> Iodine-125 Seed



Best<sup>®</sup> Palladium-103 Seed



**Do You Know?** Best Medical is also the only supplier of Iridium-192 in nylon ribbons and Gold-198 seeds!

# 24 hours

We ship your order within 24 hours, 7 days a week Anywhere in the world!

Best<sup>®</sup>  
Localization  
Needles

- Best localization needles are available in 5, 7, 10, 15 and 20 cm lengths
- Available in all ranges of seed activities (0.08 to 0.4 mCi) with a shelf life of 90 days
- Needles are supplied sterile and are individually packaged and shielded

© 2014 Best Medical International, Inc.

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
703 451 2378 800 336 4970 [www.bestmedical.com](http://www.bestmedical.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA



# Best® Total Solutions™ for Prostate Cancer Diagnosis & Treatment

## Best® Sonalis™ Ultrasound Imaging System

The Best® Sonalis™ Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimulView™ Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.

- Totally sealed, self-healing antimicrobial keyboard with SensoFoil® Technology
- Longitudinal array provides for 140 mm length of view, encompassing the bladder, prostate and perineum
- Simultaneous imaging of transverse and sagittal planes
- PC Based System provides a platform for future upgrades and application-specific modules

### Advanced Keyboard Technology



Best Medical International has been serving the Radiation Therapy Community since 1977!



## Best® UCT-1000 Universal Chair/Table

- 4-Function Hand-Held Pendant (Control)
- Electronic variable height control/back support/foot support & horizontal tilt 12 degrees – head or foot
- Padded variable height intravenous arm rest
- Removable Stainless Steel IV Pole
- Stainless steel bilateral cushioned knee supports

**Best nomos®**

Best NOMOS One Best Drive, Pittsburgh, PA 15202 USA  
412 312 6700 800 70 NOMOS [www.nomos.com](http://www.nomos.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA




## Best® Stepper/Stabilizer

- Patented cradle design provides for the absolute calibration of the transverse ultrasound image
- Infinite travel in all six axes allows for smooth, one-handed positioning without requiring fine adjustments
- Fully autoclavable for rapid OR turnaround
- Mounts to multiple-width tables and all standard rail dimensions

## BATCAM™ Multiprobe Image-Guided Radiation Therapy

BATCAM™ uses the latest compact ultrasound technology to provide increased sensitivity and resolution which yields more clinical information and expanded applications. Plus, it provides the ability to orient the probe in any rotation or translation position.

- Imports treatment plans from all major treatment planning systems
- Enhanced positioning precision using triangulated positioning techniques
- ImageSync™ Interactive Alignment provides improved visualization
- Position verification with image capture ensures accuracy
- User-friendly software components featuring a touch-screen interface



## Best® Treatment Planning System (TPS)

Best® TPS is the planning system of choice for enhanced patient outcomes.

- Dual Activity Module creates treatment plans using multiple activity sources
- Octant Therapy™ Module allows view of DVH values of specific prostate volumes
- Slice Shifter™ provides ability to correct image slice positioning errors
- Concurrent 2D & 3D Visualization allows instant visualization of seed placement, dose distribution & anatomical structures
- Pattern Load Feature – User-definable seed patterns can be used to initiate new plans

© 2014 Best Medical International, Inc.

Specifications within are subject to change.





# Best® TargetScan Touch™ — An Ultrasound-Based Stereotactic Biopsy Guidance and Tumor Localization System

## Best® Total Solutions™ for Prostate Cancer Diagnosis & Treatment

The Best® TargetScan Touch™ provides active surveillance of slow growing tumors, as well as predictive modeling for the treatment planning of more aggressive cancers.

- Point and click interface: 'Touch' the needle and instantly image in the other plane
- No need to purchase a \$20,000 digital stepper — connect straight into your treatment planning software!
- Real-time dosimetry version with two monitor mounts for system and laptop
- Smaller footprint for your Brachytherapy needs



## Best® Fiducial Marker Kit



Best Medical can provide any standard or customized marker configuration made to your exact specifications, sterile or non-sterile, shipped within 24 hours!

[www.BestTotalSolutions.com](http://www.BestTotalSolutions.com)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
703 451 2378 800 336 4970 [www.bestmedical.com](http://www.bestmedical.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

