




Best™ medical international

Total Solutions™ for Breast Cancer Treatment & Diagnosis from TeamBest Companies



Best™ Breast Localization Needles



- Best localization needles are available in 5, 7, 10 and 15 cm lengths
- Available in various ranges of seed activities with a shelf life of 90 days
- Needles are supplied sterile and are individually packaged and shielded

Best™ HDR Afterloader

- Provides 18 channels for dose delivery
- Battery back-up in case of power failure and automated wire recovery
- Source wire tracking using sensors and encoders
- Treatment data maintained in case of failure and treatment resumed after fix
- Quick source replacement process thereby reducing down time
- Available with Cobalt-60 or Iridium-192 sources. Also, Ytterbium-169 source available soon!



Best™ Double-Balloon Breast Brachytherapy Applicator



- Improved dose distribution and conformity
- Less dose to critical organs such as skin, lung, heart, chest wall, etc.
- Convenient to use
- US FDA 510(k) registered

Certain products shown may not be available for sale currently.

Best Medical International, Inc. © 2015–2018

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA