

Best[®] medical canada

MOSFET Patient Dose Verification Systems

Technology Built on 25 Years of Dosimetry Experience

▼ Portable DOSIMETER

- Stand-alone dosimetry system
- Portable & compact with LCD display
- Optional data export via USB
- No software or computer necessary



▲ mobileMOSFET Wireless

- Dose obtained in Real-Time
- Simplifies dosimetry & minimizes QA time
- Wireless with Bluetooth technology
- Export dosimetry data to R&V system

**IMRT • IGRT/TOMOTHERAPY • VMAT •
TBI • RADIOLOGY • BRACHYTHERAPY**

► MOSFET Dosimeter Features:

- Isotropic ($\pm 2\%$ for 360°)
- Active area of 0.2×0.2 mm
- Pinpoint measurement for high gradient fields
- Same calibration factor for photons & electrons at MeV energies



© 2012 Best Medical International, Inc.

Best[®] medical canada 413 March Road, Ottawa, Ontario K2K 0E4 Canada
phone 613 591 2100 877 668 6636 fax 613 596 5243 www.bestmedical.ca

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

