


# Best™ *ABT Molecular Imaging*

Best 7.5 MeV Cyclotron BG-75 BioMarker Generator with on-board chemistry & QC for in-house production of <sup>18</sup>F-FDG & other PET agents

TEAMBEST GLOBAL



[www.teambest.com](http://www.teambest.com)



The BG-75 Biomarker Generator integrates a compact mini-cyclotron, kit based micro-chemistry, and automated quality control, simplifying in-house production of <sup>18</sup>F-FDG and advanced biomarkers.

- Push button graphic interface
- Kit based chemistry
- Single or batch dose production
- Final dose delivery to syringe or vial (option)
- Automated quality control testing
- Integrated cyclotron & chemistry self-shielding
- Complete production lab in a 5 x 5 meter area

TeamBest Companies © 2020

TeamBest ASIA Private Limited email: [asia@teambest.com](mailto:asia@teambest.com)  
2/299A Plot No. 36, 5th Cross Street • V.G.P. Prabhu Nagar, Perumbakkam • Chennai – 600100 India  
[www.teambestasia.in](http://www.teambestasia.in) [www.bestabt.com](http://www.bestabt.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA