

**New!**

# A Breakthrough in Ultrasound Technology – Best® Sonalis™ Imaging System


## Featuring SimulView™ Ultrasound Technology

The Best® Sonalis™ Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimulView™ Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.

- ▶ Totally sealed, self-healing antibacterial keyboard with SensiFoil™ Technology
- ▶ Longitudinal array provides for 140 mm length of view, encompassing the bladder, prostate and perineum
- ▶ Simultaneous imaging of transverse and sagittal planes
- ▶ PC Based System provides a platform for future upgrades and application-specific modules



### Advanced Keyboard Technology



© 2010 Best Medical International, Inc.

**Best nomos** One Best Drive, Pittsburgh, PA 15202 USA  
phone 412 312 6700 800 70 NOMOS fax 412 312 6701 www.nomos.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

**Best**

healthcare for everyone