

# Best<sup>®</sup> medical *international*

## Dummy Ribbons



**Best Medical can  
provide any standard  
or customized  
configuration made  
to your exact  
specifications, sterile  
or non-sterile,  
shipped within 24  
hours — worldwide!**

## **Contact us:**

**Best Medical International  
7643 Fullerton Road  
Springfield, Virginia 22153 USA**

**Tel: 703-451-2378 • 800-336-4970**

**Fax: 703-451-5228**

**[www.bestmedical.com](http://www.bestmedical.com) • [www.teambest.com](http://www.teambest.com)**



## Mission Statement

*Everyone deserves the best health care. Our goal is to work with medical professionals to provide the Best products and services. As a leader in innovative technology, our mission is to uphold our reputation for excellence in the health care field by developing, manufacturing, and delivering cost-effective, high quality products to benefit patients around the world.*

## Corporate Values

*Fulfilling our mission requires a dedicated employee team—one that is committed to uncompromised quality and safety, and fostering a sense of pride in the company.*

## Best employees will:

- *Honor the dignity and uniqueness of each person, including fellow employees, customers, distributors, and vendors, through their actions, words, and behavior. They will treat others with respect, as they would wish to be treated.*
- *Embrace the spirit of cooperation and teamwork for which Best Medical is known, an environment which supports and recognizes effort and results.*
- *Promote and practice open and honest communication at all times.*
- *Listen to others and remain open to suggestions at all times, as we focus on achieving excellence in our work.*
- *Attract and retain leaders with exceptional abilities and integrity.*
- *Ensure all resources are managed in a responsible and effective manner.*

# Best ADVANTAGE

- Best Medical International is the ONLY supplier of Iridium-192 and Gold-198 seeds.
- We are the first and only provider of Iridium-192 in nylon ribbons in a wide range of activities.
- We are the first and only company to provide custom service for Gold-198 permanent implants, as well as a wide range of activities.
- We are the first and only company that will call our customers the morning after shipping to inform them of mode of transportation and expected arrival time.
- We can supply Ir-192, Au-198, Pd-103, I-125 and Brachytherapy accessories on a single Purchase Order.
- We ship most materials within 24 hours of receipt of the order.

---

# Catalog 2013 **Table of Contents**

Nylon and Cable Dummy Ribbons ..... 1





## Anderson Forell Suthanthiran Nylon Dummy Ribbon Set

Part # 5006-0003

### Set Contains:

One Ribbon with 12 Seeds each of 1.2 mm Length

One Ribbon with 12 Seeds each of 1.8 mm Length

One Ribbon with 12 Seeds each of 2.4 mm Length

One Ribbon with 12 Seeds each of 3.0 mm Length

One Ribbon with 12 Seeds each of 4.0 mm Length

This enables identification of ribbons in complex implants, facilitates implant radiography and improves localization accuracy. Dr. Lowell L. Anderson of Memorial Sloan-Kettering Cancer Center and Mr. Bruce Forell of City of Hope Medical Center have been very instrumental in bringing about this design concept.

## Nylon Dummy Ribbon

Part # 220



The diameter of the Nylon Dummy Ribbon and High Density Seed are 0.85 mm and 0.5 mm respectively. Nylon Dummy Ribbons are also available with seed lengths of 0.6 mm, 1.2 mm, 1.8 mm, 2.4 mm, 3.0 mm, 4.0 mm and 5.0 mm.


## Cable Dummy Ribbon

Part # 423/HDR Part # 530



Cable Dummy Ribbons are similar to the Nylon Dummy Ribbons. The diameters of the cable and the dummy seeds are 0.5 mm and 0.85 mm respectively. Custom made Cable Dummy Ribbons are available with or without alternate markers using seed lengths of 0.6 mm, 1.2 mm, 1.8 mm, 2.4 mm and 3.0 mm.

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.



**Best Medical International, Inc.** as author/distributor, believes the technical information and statements contained in this catalog to be reliable and correct at the time of publication. However, Best Medical International, Inc., does not assume any liability for statements, information or recommendations based on/or made in reliance on the information contained in this publication.

Data provided herein are for informational purposes only and are not intended for clinical use of any nature. Use of these data in such a manner is the sole responsibility of the user, and the author/distributor of the information makes no claim or inference as to its clinical suitability and/or accuracy. Best Medical International, Inc. assumes no warranty or liability for use of this information or products. Designs and specifications are subject to modification and/or change(s).

© 2012 Best Medical International, Inc.

## Terms and Conditions

**Ordering:** Orders may be placed with a purchase order number and transmitted by FAX, phone or e-mail.


**Payment:** Within 30 days.

**Freight:** FOB origin, freight prepaid and charged back through invoice.

**Returns:** Unused accessories may be returned to Best Medical International, Inc. within 15 days if found unsuitable for a particular application. No return or authorization number is necessary. A note stating the reason for return, name, and address of hospital will expedite processing of any credit.

**Note:** Best Medical International, Inc. assumes no liability for consequential damages of any kind for the material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer.

Best Medical International, Inc., reserves the right to modify the dimensions of accessories/components as well as the prices.



# Are you Searching for the Greatest Healthcare Company in the World?



**Look to TeamBest®!**

World population is increasing by more than 100 million per year. To meet the needs of our aging and growing population, TeamBest® companies are committed to making quality healthcare affordable and accessible globally.

**Ready to Serve Globally for All Your Needs**

- Arplay Medical
- Best Cyclotron Systems
- Best Medical Africa
- Best Medical Asia
- Best Medical Canada
- Best Medical China
- Best Medical Engineering
- Best Medical Europe
- Best Medical International
- Best Medical Latin America
- Best Medical Middle East
- Best Medical Parts
- Best Medical Research
- Best Medical Russia
- Best NDT
- Best NOMOS
- Best Particle Therapy
- Best Theratronics
- Best Vascular (Novoste)
- CNMC Company
- Genetherapy Pharmaceuticals
- Huestis Diagnostics
- Huestis Medical and more...

**For more information, visit [www.teambest.com](http://www.teambest.com)**

© 2013 Best Medical International, Inc.



Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
703 451 BEST (2378) 800 336 4970 [www.bestmedical.com](http://www.bestmedical.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA