

Best[™] nomos[®] Cyber Sonalis[®] Ultrasound System

The Best[™] Cyber Sonalis[®] Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimulView[™] Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.



Additional Features

- ✓ Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- ✓ 21.5" touch screen technology
- ✓ Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- ✓ Patented probe design
- ✓ Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- ✓ Advanced modular software design provides for future upgrade path
- ✓ Convertible to desktop or portable without losing any functionality
- ✓ Dual battery backup built-in to protect from power loss



**Certain products shown are not available for sale currently.*

Sonalis[®] Transducers

We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the **BEST** name in ultrasounds!

8L2A Linear Array



15L4A Linear Array



9MC3 Curved Array



8BP4 Bi-Plane



12L5A Linear Array



16L5 Linear Array



8EC4A Endocavity



8TE3 Trans-esophageal



14L3 Linear Array



8V3 Phased Array



XY-BI-Plane Phased Array



Pedoff



15LW4 Linear Array



4V2A Phased Array



10EC4A Endocavity



16HL7 High Frequency Linear Array



15LA Linear Array



5C2A Curved Array




10BP4 Bi-Plane



Best Nomos • One Best Drive • Pittsburgh, Pennsylvania 15202 USA
tel: 412 312 6700 800 70 NOMOS www.nomos.com

Best™ Integrated Radiation Therapy Solutions




Best™ Flexi Needles



Best™ Stepper/Stabilizer



Best™ Brachytherapy Kit



* Half Moon Buttons are available in either radiopaque (pictured above) or clear nylon.
**Catheter tubes are available in purple, green, yellow, clear or blue with either radiopaque or clear nylon buttons.


BDS Whole Body Dosimeters



BDS Extremity Ring Dosimeters



Best™ X-ray Marker Cables



Best™ ABT Molecular Imaging



Best™ B35ADP

Alpha/Deuteron/Proton Cyclotron

- **Proton Particle Beam:**
1000 μ A Beam Current
up to 35 MeV Energy
- **Deuteron Particle Beam:**
500 μ A Beam Current up
to 15 MeV Energy
- **Alpha Particle Beam:**
200 μ A Beam Current
up to 35 MeV Energy